

Sample Fixation Protocol for Pico/Nanoplankton Flow Cytometry

Preparation of Paraformaldehyde Solution

Materials

Paraformaldehyde powder
1N NaOH
stirring/hot plate
Chemical fume hood
pH meter

Mix 900 mL DI water and 100 g paraformaldehyde powder. Set up on a stirring/hot plate under hood. Heat to approximately 60 °C. Do not boil! Stir for approx 1hr. Turn off heat. Add 100 μ L 1N NaOH to “clear” solution. Cool to room temperature. Add 100 mL phosphate buffer solution or filtered seawater. Filter through GF/F filter to remove precipitate. Test pH. Should be 7.4 - 8.0 (approx. equal to seawater). If necessary, add more NaOH. This yields a 10% solution (approximately).

Sample Fixation Procedure

Materials

cryovials with O-rings (eg., Wheaton #985747 Cryule vials)
liquid N₂
aluminum canes
GF/F filters and filtration apparatus
pipettors

1. Add 50 μ L of 10% paraformaldehyde to a labeled cryovial.
2. Add 1 mL sample to cryovial, cap, and mix. This yields a 0.5% final conc.
3. Allow fixation for 1 – 8 h in refrigerator (4°C).
4. Place vials in labeled aluminum canes and put directly into liquid N₂.
5. Store samples in liquid N₂ – can be transferred to –80°C freezer.

Notes:

1. If sample has large particles (detritus, zooplankton) pre-screen through 100 – 150 μ m mesh nytex.
2. Cryovials can be pre-loaded with paraform., then add sample, etc.
3. Method is derived from Vaultot et al. 1989 (Vaultot, D., Courties, C. & Partensky, F. 1989. A simple method to preserve oceanic phytoplankton for flow cytometric analyses. *Cytometry* 10:629-635).